



## TECHNICAL DATA

### Nacera® Hybrid

Type: Blanc, Block M, Block S

Color: A1, A2, A3, B1, C2, D2, BL2, BL3.5

### PHYSICAL PROPERTIES

Property	Values
Flexural Strength	175 MPa
Compressive Strength	490 MPa
Modulus of Elasticity	9900 MPa
Vickers Hardness HV 10	700 MPa

Notified Body: ecm GmbH **CE 0481**

These properties were measured on test samples. The values are typical material properties and may vary according to product configuration, geometry and manufacturing process.

2017-09-06